


Filter set 02
488002 - 0000

(typical curves)

— excitation: G 365
— beamsplitter: FT 395
— emission: LP 420



**Filterset 02 for the following
Fluorochromes:**

- 1,5 IAEDANS
- 1,8-ANS
- 4-Methylumbelliferone
- 7-Amino-4-methylcoumarin
- 7-Hydroxy-4-methylcoumarin
- ABQ
- Alexa Fluor 350™
- AMC, AMCA-S
- AMCA (Aminomethylcoumarin)
- AMCA-X
- Aminocoumarin
- Aminomethylcoumarin (AMCA)
- Anthrocyll stearate
- Bisbenzamide
- Bisbenzimidide (Hoechst)
- Calcein Blue
- Cascade Blue™
- Dansyl

→ Fluorochromes

📄 pdf Filterset

Dansyl Amine
Dansyl Cadaverine
Dansyl Chloride
Dansyl DHPE
DAPI
Dapoxyl 2
Dapoxyl 3
DIDS
Dopamine
ELF 97
Fast Blue
FluoroGold (Hydroxystilbamidine)
Granular Blue
Hoechst 33258
Hoechst 33342
HPTS
Hydroxycoumarin
Hydroxystilbamidine (FluoroGold)
Indo-1, low calcium
Intrawite Cf
Laurodan
Leucophor PAF
Leucophor SF
Lyso Tracker Blue
LysoSensor Blue
LysoSensor Yellow/Blue
Marina Blue
Noradrenaline
Nuclear Yellow
PhotoResist
Phorwite AR
Phorwite BKL
Phorwite Rev
Phorwite RPA
PhotoResist
Pyrozal Brilliant Flavin 7GF
Serotonin
sgBFP™ (super glow BFP)
SITS
SITS (Stilbene Isothiosulphonic Acid)

SPQ (6-methoxy-N-(3-sulfopropyl)

quinolinium)

Stilbene

Thioflavin TCN

Thiolyte

Tinopol CBS (Calcofluor White)

True Blue

Uvitex SFC

Y55F

Y66H