


# Filter set 10

488010 - 0000

(typical curves)

— excitation: BP 450 - 490  
 — beamsplitter: FT 510  
 — emission: BP 515 - 565



Filterset 10 for the following Fluorochromes:

- 5-Carboxyfluorescein (5-FAM)
- 5-FAM (5-Carboxyfluorescein)
- Acridine Orange, both DNA & RNA
- Acridine Yellow
- Alexa Fluor 488™
- ATTO-TAG™ CBQCA
- Auramine
- Aurophosphine
- BCECF (high pH)
- BCECF (low pH)
- Calcein
- CFDA
- CMFDA
- CyQuant
- DTAF
- FDA
- FITC
- FITC Antibody
- Fluo-3

→ Fluorochromes

PDF Filterset

Fluo-4
Fluorescein (FITC)
Fluorescein Diacetate
Fluoro-Emerald
FluorX
GFP (S65T)
GFP red shifted (rsGFP)
Lyso Tracker Blue-White
Lyso Tracker Red
Mitotracker Green FM
Oregon Green 488-X
Oregon Green™ 488
Oregon Green™ 500